Description:

The LaserDock is an accessory to be used with surgical lasers and laser fibers. The device is a sterile temporary holding device for fiber optic laser energy delivery devices (laser fibers) when the laser fiber is not in use during the clinical procedure. The LaserDock securely holds the laser fiber in place during clinical procedures when the fiber is not deployed in the patient. This is a single use device which consists of an inner medical grade silicone tube for resting the laser fiber inside and preventing any inadvertently discharged laser energy from escaping and injuring the users and patients, a clamp for securing the laser fiber in place, and clam shell housing to protect the inner silicone tubing.

Indications for Use:

The LaserDock is designed for use with Ho:YAG, Nd:YAG, CO2 and KTP laser systems when laser fibers are deployed.

Warnings:

- **Caution: Federal law restricts this device to sale by or on the order of a physician.**
- Always consult and follow the laser or laser fiber manufacturer's package insert or IFU and read all warnings and cautions before using a laser or laser fiber equipment. Read all directions for use to ensure safe and proper operation.
- Should the laser energy accidentally be discharged while the laser fiber is inside the LaserDock, the device must be discarded and a new LaserDock used.
- Should the laser energy accidentally be discharged while the laser fiber is inside the LaserDock, the laser fiber must be discarded and a new laser fiber used.
- Should the laser fiber break while in the LaserDock, a new laser fiber should be opened for use in the patient, or the end of the laser fiber can be "stripped" according to standard technique to prepare the laser fiber for re-use in the patient.
- Should the laser fiber break while in the LaserDock, the LaserDock should be discarded and a new one opened for use.
- Accidental laser energy discharge while the laser fiber is inside the LaserDock may produce a loud noise and visible smoke. Should this occur, replace both the laser fiber and the LaserDock. The LaserDock is designed to protect the user and patient from injury due to accidental discharge of the laser fiber while housed in the LaserDock. There is no risk to the user or patient if this occurs.

Caution: Laser fiber should be inserted approximately 8cm (3 in) into the LaserDock (see Figure 1). An arrow and line mark on the side of the LaserDock outer housing indicate this minimum distance. Locate the distal end of the laser fiber against the indicator mark and measure back to the opening of the LaserDock to determine the minimum length of fiber to be inserted. This will ensure the laser fiber is inserted sufficiently into the silicone tubing in case the laser fiber accidentally discharges while in the device, so that any laser energy discharge will be properly contained inside the device. Do not over-tighten the LaserDock and therefore damaging the largest fiber sizes, or under-tighten the LaserDock providing insufficient retention strength on the smallest laser fiber sizes.

The LaserDock is a single-use device and may not be re-used or re-sterilized.

Precautions:

- Use the LaserDock only with the appropriate Personal Protective Equipment (PPE) as recommended by the laser or laser fiber manufacturer.
- Use of the LaserDock does not substitute or eliminate the need to follow all laser fiber safety precautions. The LaserDock is designed to add safety when handling laser fibers in a clinical setting, but is not meant to supersede any warnings or precautions contained in the package insert or IFU of the laser and/or laser fiber device being used.
- Use caution when using laser fibers in high oxygen concentration environments where there is a risk of igniting flammable accessories or equipment (bronchoscopes, surgical drapes absorbent packing, etc.);
- Be careful to not allow stray radiation or reflected laser beams/energies during surgery or clinical procedures;
- There are risks to tissues or organs which may be adversely affected if struck by stray or scattered laser radiation during surgery; and
- Steps should always be taken to protect patients from stray laser radiation.
- LaserDock was validated to be used with laser fibers up to 2.3mm in diameter.

Directions for Use:

1. Open the LaserDock packaging using sterile technique and remove the device.
2. Place the flat, raised side of the LaserDock down on the surgical field
3. Secure the LaserDock on the surgical field by clamping a surgical clamp to a surgical drape onto each of the LaserDock "handles". These appear like "protrusions" on the external housing of the device and are used to assist in securing the device to the surgical drape. The surgical clamp should face away from the LaserDock latch clamp so as to avoid interference with laser fiber insertion into the LaserDock. See Figure 1.
4. Gently rotate the LaserDock latch clamp to the right and to the left to verify correct latch operation. Clockwise rotation tightens the clamp whereas counter-clockwise will loosen it (see Figure 2). Please see following section for detailed instructions on locking and unlocking the laser fiber wire in the clamp.
5. Keep the clamp in the open position until the laser fiber is fully inserted into the LaserDock.
6. Open the laser fiber package according to manufacturer's instructions and verify its integrity before inserting it into the LaserDock.
7. Place the operative end of the laser fiber into the opening of the LaserDock and manually advance it gently approximately 8cm (3 in) past indicator mark until resistance is met. Note that there is a bend in the tubing which indicates to the user that the laser fiber has been fully inserted. See Figure 2.
8. Turn the latch clamp slowly to the right, until resistance is met, to lock the laser fiber in place.
9. To release the laser fiber for use, rotate the latch clamp gently to the left and carefully remove the laser fiber from the LaserDock.
10. The laser fiber can be inserted and the latch clamp rotated to lock the laser fiber in the LaserDock as often as needed during the surgical case.
11. The laser fiber should always be inspected for integrity before each use in the patient.
12. The laser energy source should always be moved to the "standby" position when the laser fiber is not in active use in the patient and when inside the LaserDock.

Figure 1



Figure 2

Manufactured by:
Patient Pocket, LLC
405 Hialeah Drive
Cherry Hill, NJ 08002
Phone: (216) 272-0805

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Medes Limited